	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.0
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2022.11.01

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Transpro CD 01 is a universal transient medium, which can be used for subculture, high-density culture and transient transfection culture of HEK293 cells and CHO cells. The transient transfection process does not require centrifugation to change the medium. Transpro CD 01 is suitable for the use of 293 series cells such as HEK 293, Expi293F, 293F, 293E and CHO series cells such as ExpiCHOS and CHOS for transient transfection expression culture of antibodies, recombinant proteins and viruses during the development and manufacture process. Transpro CD 01 is an animal-derived component free (ACF), protein free (PF), chemically defined (CD) medium. Transpro CD 01 medium does not contain any growth factor and hydrolysates, which ensures consistency between batches and improves the efficiency of the cell culture process. Transpro CD 01 liquid package does not contain L-glutamine, HT and anti-clumping agent.

### 【细胞培养 Cell culture】

- ① 建议细胞接种密度 Suggested cell inoculation density:  $0.2\sim 1.0\times 10^6$  cells/mL
- ② 温度 Temperature: 37°C
- ③ CO<sub>2</sub>: 8%
- ④ 使用时需要额外添加 L-谷氨酰胺 4-6mM

It needs to be supplemented with 4-6mM L-glutamine when used.

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Most cells lines can adapt directly into this product. They can be directly inoculated into this medium and passed more

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.0
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2022.11.01

than three times. For some cell lines, sequential cell adaptation may be used when using this medium.

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Prepare the cryopreservation solution on the super clean workbench: 90% Transpro CD 01 + 10% dimethyl sulfoxide (DMSO) mixture, precooling at 2~ 8°C(Temperature will be released when DMSO is diluted);

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Cryopreserved cell fluid: in logarithmic growth stage, with a density greater than  $1.5 \times 10^6$  cells / ml, and the activity rate is greater than 95%. Generally, the recommended freezing density is  $1.0 \sim 1.5 \times 10^7$  cells/mL;

- ③ 细胞液 800rpm 离心 5 min； The cell fluid was centrifuged at 800 rpm for 5 min;

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Slowly pour out the supernatant and resuspend the cells with cryopreservation solution, and the cryopreservation density is  $1.0 \sim 1.5 \times 10^7$  cells / ml, transfer the cells to the sterile cryopreservation tube;

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at - 80 °C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4℃冻存 30min；
- freeze at 4°C for 30min；
- -20℃冻存 2~4 小时；
- freeze at -20°C for 2~4h；
- -80℃冻存过夜；
- freeze at - 80°C overnight；
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare a 37 °C water bath to thaw frozen cells ;

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15 ml sterile centrifuge tube is prepared, and 2~5mL Transpro CD 01 is added;

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Take out the cryopreservation tube from the liquid nitrogen tank and rapidly thaw (<1 minute) frozen cells in a 37°C water bath;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell suspension to a 15 ml centrifuge tube and centrifuge at 800 rpm for 5 min;

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Slowly pour out the supernatant, resuspend it with 15 ~ 20 ml preheated Transpro CD 01, and transfer it to a 125 ml shake flask;

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Place it in a shaking incubator with 8% CO<sub>2</sub>，110 ~ 130rpm，at 37°C for culture;

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After 2-3 days of culture, the cells are counted and subcultured.

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The cells are seeded at 0.2 ~ 1.0 ×10<sup>6</sup> cells/ml, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. When the cell viability returned to normal and reached more than 90%, it was increased by 0.2 ~ 1.0×10<sup>6</sup> cells/ml was amplified until the required seed volume and normal seed state were reached: the vitality was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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The day before transfection need to seed cells at 2.0×10<sup>6</sup> viable cells/ml, the cell density can reach 4.0×10<sup>6</sup> viable cells/ml on the second day ;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.0
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2022.11.01

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After cell counting on the second day of culture, the cell viability was more than 95%, and the living cell density was  $\geq 4.0 \times 10^6$  cells / ml, can be used directly; If the cell density is lower than  $4.0 \times 10^6$  cells / ml, the cells can be collected by centrifugation (800 rpm, 5 min), and the cells can be separated at  $4.0 \times 10^6$  cells / ml density was resuspended in transpro CD01 medium;

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The mixture of DNA and PEI was prepared according to the optimized transient process;

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Add the mixed solution to the culture medium for culture;

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After 18 hours of culture, it is recommended to supplement the supplemented medium Transpro feed 1 (the concentration is recommended to be 3~5% of the initial culture volume), or the combined supplemented medium DN feed B2 (the concentration is recommended to be 0.3~0.5% of the initial culture volume), which can further improve the density of living cells and protein expression.

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Culture until 7 days, or the vitality is less than 60%, and end the culture.

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Liquid packaging: 2°C to 8°C, Protect from light; Shelf life: 12 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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网址 Website: www.duoningbio.com